

CHOLESTEROL

Enzymatic-spectrophotometric
CHOLESTEROL OXIDASE/PEROXIDASE

The application parameters comprised here constitute a guide to facilitate the validation of our reagents by the instrument. It is advisable to validate the use when there is any change in software or reagent versions.

Instrument: ARCHITECT C8000

Reagent preparation

Reagent is ready to use.

Instrument settings

General				Conventional Units		
Assay Number	CHOL ...	Type Version	Photometric 1	Results Unit	mg/dL	
				Assay Defaults	Low Linearity	0.3
					High Linearity	1000
Reaction Definition				<u>Gender and Age specific Ranges</u>		
Reaction Mode	End Up		Read Times	GENDER	AGE (Y)	NORMAL
Wavelength	500 / 0	Main	16 -- 18	Male	0-130	0 - 239
Last Required Read	18	Flex	— -- —	Female	0-130	0 - 239
Absorbance Range	... / ...	Color Correction	— -- —	Either	0-130	0 - 239
Sample Blank Type	None	Blank Read Time	— -- —			
Reagent / Sample				Results Units	mg/dL	
Reagent	CHOL	Reagent Volume	R1 R2 300	Decimal Places	1	(Range 0-4)
Diluent	Saline	Water Volume	— -- —	Correlation Factor	1.0000	
Diluent Disp. Mode	Type 0	Dispense Mode	Type 0 Type 0	Intercept	0.0000	
Diluent Name	Sample	Water	Dilution Factor			
Standard	3		1:1.00			
Validity Checks						
Reaction definition	None					
Rate Linearity %	—					
Calibration						
Method	Linear					
Calibrators and Volumes						
Calibrator Set	Standard			Calibrator Level	Concentration	Volume
				Blank	0.0	3
				Std. 1	(*)	3
				Std. 2	--	--
				Std. 3	--	--
				Std. 4	--	--
				Std. 5	--	--
Replicates	3	(Range 1-3)				
Intervals				Validity Checks		
Calibration Interval	Full Interval	999 hours		Blank Absorbance range:	0.0000 -- 0.0000	
Calibration Type	Adjust Type	None		Span Absorbance Range:	0.0000 -- 0.0000	
				Expected cal Factor :	0.00	
				Expected Cal factor tolerance:	% 0	
(*) Enter the standard value						
... Values to be entered by the user						