

TRANSFERRIN

Turbidimetry

The application parameters comprised here constitute a guide to facilitate the validation of our reagents by the instrument. It is advisable to validate the use when there is any change in software or reagent versions.

Instrument: ARCHITECT C8000

Reagent preparation

Reagent is ready to use.

Instrument settings

General				Conventional Units		
Assay Number	TF ...	Type Version	Photometric 1	Results Unit	mg/dL	
Reaction Definition				Assay Defaults	Low Linearity	4.8
Reaction Mode	End Up		Read Times	High Linearity	700	
Wavelength	340 / 700	Main	14 -- 16	Gender and Age specific Ranges		
Last Required Read	16	Flex	— -- —	GENDER	AGE (Y)	NORMAL
Absorbance Range	... / ...	Color Correction	— -- —	Male	0-130	200 - 360
Sample Blank Type	None			Female	0-130	200 - 360
Reagent / Sample				Reults Units	mg/dL	
Reagent	TF	Reagent Volume	R1 R2 200	Decimal Places	1	(Range 0-4)
Diluent	Saline	Water Volume	— —	Correlation Factor	1.0000	
Diluent Disp. Mode	Type 0	Dispense Mode	Type 0 Type 0	Intercept	0.0000	
Diluent Name	Sample	Water	Dilution Factor			
Standard	2		1:1.00			
Validity Checks						
Reaction definition	None					
Rate Linearity %	—					
Calibration						
Method	Logit - 4					
Calibrators and Volumes						
Calibrator Set	Standard		Calibrator Level	Concentration	Volume	
			Blank	0.0	2	
			Std. 1	(*)	2	
			Std. 2	(*)	2	
			Std. 3	(*)	2	
			Std. 4	(*)	2	
			Std. 5	(*)	2	
Replicates	3	(Range 1-3)				
Intervals				Validity Checks		
Calibration Interval	Full Interval	999 hours	Blank Absorbance range:	0.0000 -- 0.0000		
Calibration Type	Adjust Type	None	Span Absorbance Range:	0.0000 -- 0.0000		
				Expected cal Factor :	0.00	
				Expected Cal factor tolerance:	% 0	
(*) Enter the standard value						
... Values to be entered by the user						