

α-AMYLASE

Continuous-spectrophotometric
IFCC

The application parameters comprised here constitute a guide to facilitate the validation of our reagents by the instrument. It is advisable to validate the use when there is any change in software or reagent versions.

Instrument: AU 400 / 640

Reagent preparation

Working Reagent: Pour the contents of the Reagent B into the Reagent A bottle. Mix gently. Other volumes can be prepared in the proportion: 4 mL Reagent A + 1 mL Reagent B. Stable for 20 days at 2-8°C.

Instrument settings

Specific Test Parameters		Min OD	Max OD
Sample: Volume	6 µL	Dilution	0 µL
Reagents: R1 Volume	200 µL	Dilution	0 µL
Wavelength: Pri.	410	Sec.	480
Method:	RATE		
Reaction slope:	+		
Measuring Point 1: First	3	Last	12
Measuring Point 2: First	0	Last	0
Linearity:	15 %		
No Lag Time:	NO		
		L 0.0	H 1.5
		Reagent OD limit:	
		First L -0.1	First H 1.2
		Last L -0.1	Last H 1.2
		Dynamic Range:	
		L 0	H 1300
		Correlation Factor:	
		A 1	B 0
		On-board stability period:	
7. None selected		L	H
8. Out of range		28	100
Panic Value	L	H	
		Unit: U/L	Decimal places: 0
Calibration Specific			
Calibration Type:	MB	Formula: Y=AX+B	Counts: 2
			Process: CONC.
Point 1:	Cal. N° #	OD	CONC *
			Factor/OD-L -9999
			Factor/OD-H 9999
MB Type Factor:	5398		
* Enter the assigned standard value			
# User defined			