

# PROTEIN

Spectrophotometric  
BIURET

The application parameters comprised here constitute a guide to facilitate the validation of our reagents by the instrument. It is advisable to validate the use when there is any change in software or reagent versions.

## Instrument: **AU 400 / 640**

### Reagent preparation

Reagent is ready to be used.

### Instrument settings

Specific Test Parameters					
Sample:	Volume	4 $\mu$ L	Dilution	0 $\mu$ L	
Reagents: R1	Volume	200 $\mu$ L	Dilution	0 $\mu$ L	
Wavelength:	Pri.	540	Sec.		
Method:	END				
Reaction slope:	+				
Measuring Point 1:	First	0	Last	27	
Measuring Point 2:	First	0	Last	0	
Linearity:	%				
No Lag Time:					
			Reagent OD limit:	First L	-0.5
				First H	0.1
				Last L	-0.5
				Last H	0.1
			Dynamic Range:	L	0
				H	150
			Correlation Factor:	A	1
				B	0
			On-board stability period:		
			L	H	
7. None selected					
8. Out of range			65	80	
Panic Value			L	H	
			Unit: g/L	Decimal places: 0	
Calibration Specific					
Calibration Type:	AB	Formula:	Y=AX+B	Counts:	2
				Process:	CONC.
Point 1:	Cal. N°	OD	CONC	Factor/OD-L	Factor/OD-H
	#		*	-9999	9999
* Enter the assigned standard value					
# User defined					