

# RHEUMATOID FACTORS

Turbidimetry  
LATEX

## Non linear calibration

The application parameters comprised here constitute a guide to facilitate the validation of our reagents by the instrument. It is advisable to validate the use when there is any change in software or reagent versions.

Instruments: **AU 400 / 640**

## Reagent preparation

Reagent 1: use the Diluent. Stable for 28 days at 2-8 °C.

Reagent 2: use the Latex. Shake the Latex vial gently before using.

## Instrument settings

Specific Test Parameters					
Sample:	Volume	3 $\mu$ L	Dilution	0 $\mu$ L	
Reagents: R1	Volume	240 $\mu$ L	Dilution	0 $\mu$ L	Reagent OD limit:
	R2	Volume	60 $\mu$ L		First L -2.0 First H 2.5
Wavelength:	Pri.	660	Sec.		Last L -2.0 Last H 2.5
Method:	END				
Reaction slope:	+		Dynamic Range:		
Measuring Point 1:	First	0	Last	20	L 0 H #
Measuring Point 2:	First	0	Last	12	Correlation Factor:
Linearity:	%		A 1 B 0		
No Lag Time:	On-board stability period:				
7. None selected			L	H	
8. Out of range			0	30	
Panic Value	L		H		Unit: IU/mL
Decimal places: 1					
Calibration Specific					
Calibration Type:	5AB	Formula:	POLYGONAL	Counts:	2
				Process:	CONC.
	Cal. N <sup>o</sup>	OD	CONC	Factor/OD-L	Factor/OD-H
Point 1:	#		*	-2.0	2.5
Point 2:	#		*	-2.0	2.5
Point 3:	#		*	-2.0	2.5
Point 4:	#		*	-2.0	2.5
Point 5:	#		*	-2.0	2.5
* Enter the assigned standard value					
# User defined					