

ALBUMIN

Spectrophotometric
BROMOCRESOL GREEN

The application parameters comprised here constitute a guide to facilitate the validation of our reagents by the instrument. It is advisable to validate the use when there is any change in software or reagent versions.

Instrument: AU 5200

Reagent preparation

Reagent is ready to use.

Instrument settings

Operation	YES	Normal Value	H	50.0
Sample Volume	3 μ L		L	35.0
Sample Volume Repeat Run	2 μ L	Dynamic Range	H	70.0
Reagent Volume R1	300 μ L		L	0.0
R2	0 μ L	Panic Value	H	60.0
W3 Parameters Prime (Bottle to Valve) R1	YES		L	0.0
R2		Reagent OD Range		
W3 Parameters Prime (Bottle to Valve) R1	100	First Point	H	0.500
R2			L	-0.100
Method	END	Last Point	H	0.500
Wavelength 1	600		L	-0.100
2	800	Decision Level for Repeat Run	H	60.0
Reaction Slope	+		L	0.0
Measuring Point S	0	R. Dispense by Black Rack		1
E	9	R. Dispense by Sample Interval		0
OD Value Range Max		W6 Bottle Side Cleaning		Yes
Min				
% Variance				
Calibration Type	AB			
ODO Conc.	0 g/L			
Conc.	*	* Enter the assigned standard value		