

DIRECT BILIRUBIN

Spectrophotometric
DIAZOTIZED SULFANILIC

The application parameters comprised here constitute a guide to facilitate the validation of our reagents by the instrument. It is advisable to validate the use when there is any change in software or reagent versions.

Instrument: **AU 5200**

Reagent preparation

Working Reagent: Transfer the contents of one Reagent B-D vial into a Reagent A-D bottle. Mix thoroughly.
Stable for 20 days at 2-8 °C.

Instrument settings

Operation		YES	Normal Value	H	0.25
Sample Volume		10 µL		L	0.00
Sample Volume Repeat Run		5 µL	Dynamic Range	H	15.00
Reagent Volume	R1	250 µL		L	0
	R2	0 µL	Panic Value	H	15.00
W3 Parameters Prime (Bottle to Valve)	R1	YES		L	0
	R2		Reagent OD Range		
W3 Parameters Prime (Bottle to Valve)	R1	100	First Point	H	1.000
	R2			L	-0.100
Method		END	Last Point	H	1.000
Wavelength	1	540		L	-0.100
	2	660	Decision Level for Repeat Run	H	15.00
Reaction Slope		+		L	0
Measuring Point	S	0	R. Dispense by Black Rack		1
	E	9	R. Dispense by Sample Interval		0
OD Value Range	Max		W6 Bottle Side Cleaning		Yes
	Min				
% Variance					
Calibration Type		AB			
ODO Conc.		0 mg/dL			
Conc.		*	* Enter the assigned standard value		