

PROTEIN

Spectrophotometric
BIURET

The application parameters comprised here constitute a guide to facilitate the validation of our reagents by the instrument. It is advisable to validate the use when there is any change in software or reagent versions.

Instrument: **AU 5200**

Reagent preparation

Reagent is ready to be used.

Instrument settings

Operation		YES	Normal Value	H	87
Sample Volume		4 µL		L	67
Sample Volume Repeat Run		2 µL	Dynamic Range	H	150
Reagent Volume	R1	350 µL		L	0
	R2	0 µL	Panic Value	H	120
W3 Parameters Prime (Bottle to Valve)	R1	YES		L	0
	R2		Reagent OD Range		
W3 Parameters Prime (Bottle to Valve)	R1	100	First Point	H	0.300
	R2			L	-0.400
Method		END	Last Point	H	0.300
Wavelength	1	540		L	-0.400
	2	660	Decision Level for Repeat Run	H	120
Reaction Slope		+		L	0
Measuring Point	S	0	R. Dispense by Black Rack		1
	E	16	R. Dispense by Sample Interval		0
OD Value Range	Max		W6 Bottle Side Cleaning		Yes
	Min				
% Variance					
Calibration Type		AB			
ODO Conc.		0 g/L			
Conc.		*	* Enter the assigned standard value		