

# α-AMYLASE

Continuous-spectrophotometric  
DIRECT SUBSTRATE

The application parameters comprised here constitute a guide to facilitate the validation of our reagents by the instrument. It is advisable to validate the use when there is any change in software or reagent versions.

## Instrument: AIRONE 200 (=AUTOHUMALYZER 900 S)

### Reagent preparation

Reagent is ready to use.

### Instrument settings

UNITS	U/L
DECIMALS	0
REACTION TYPE	Kinetic
VOLUME SERUM	8
PLASMA	8
URINE	
REAGENT 1	400
REAGENT 2	0
ABS RANGE MIN.	0.0
ABS RANGE MAX.	2.000
LINEARITY LIMIT	1300
REAGENT BLANK	No
CONTAMINATING	No
DIFERENTIAL	No
FILTER 1	405
FILTER 2	None
TIMES MIX 1	0.2
INCUBATION 1	40
LAG PHASE	20
MIX 2	0.0
INCUBATION 2	0
MEASURE	30
MEASURE TYPE	Factor
FACTOR	3292
NORMAL RANGE	0 – 60

Vers. 0207