

BILIRUBIN (TOTAL)

Spectrophotometric
DIAZOTIZED SULFANILIC

The application parameters comprised here constitute a guide to facilitate the validation of our reagents by the instrument. It is advisable to validate the use when there is any change in software or reagent versions.

Instrument: **AIRONE 200 (=AUTOHUMALYZER 900 S)**

Reagent preparation

Reagent 1: Use the Reagent A-T.

Reagent 2: Use the Reagent B-T.

Instrument settings

UNITS	mg/dL
DECIMALS	2
REACTION TYPE	End Point
VOLUME SERUM	30
PLASMA	30
URINE	
REAGENT 1	320
REAGENT 2	80
ABS RANGE MIN.	0.0
ABS RANGE MAX.	2.000
LINEARITY LIMIT	15
REAGENT BLANK	Yes
CONTAMINATING	Yes
DIFERENTIAL	Yes
FILTER 1	546
FILTER 2	None
TIMES MIX 1	0.0
INCUBATION 1	300
LAG PHASE	2
MIX 2	0.0
INCUBATION 2	0
MEASURE	1
MEASURE TYPE	Std
FACTOR	0.0
NORMAL RANGE	0.00 – 1.10

Vers. 0707