

CHOLINESTERASE

Continuous-spectrophotometric
BUTYRYLTHIOCHOLINE

The application parameters comprised here constitute a guide to facilitate the validation of our reagents by the instrument. It is advisable to validate the use when there is any change in software or reagent versions.

Instrument: **AIRONE 200 (=AUTOHUMALYZER 900 S)**

Reagent preparation

Working Reagent: Pour the contents of the Reagent B into the Reagent A bottle. Mix gently.
Stable for 3 days at 2-8°C.

Instrument settings

UNITS	U/L
DECIMALS	1
REACTION TYPE	Kinetic
VOLUME SERUM	6
PLASMA	6
URINE	
REAGENT 1	400
REAGENT 2	0
ABS RANGE MIN.	0.0
ABS RANGE MAX.	25000
LINEARITY LIMIT	25000
REAGENT BLANK	No
CONTAMINATING	No
DIFERENTIAL	No
FILTER 1	405
FILTER 2	None
TIMES MIX 1	0.2
INCUBATION 1	90
LAG PHASE	90
MIX 2	0.0
INCUBATION 2	0
MEASURE	30
MEASURE TYPE	Std.
FACTOR	0.0
NORMAL RANGE	3930 – 11500

Version 0709