

# FRUCTOSAMINE

Spectrophotometric  
NBT

The application parameters comprised here constitute a guide to facilitate the validation of our reagents by the instrument. It is advisable to validate the use when there is any change in software or reagent versions.

## Instrument: **AIRONE 200 (=AUTOHUMALYZER 900 S)**

### Reagent preparation

Reagent ready to be used.

### Instrument settings

UNITS	mmol/L
DECIMALS	2
REACTION TYPE	Kinetic
VOLUME SERUM	20
PLASMA	20
URINE	
REAGENT 1	400
REAGENT 2	0
ABS RANGE MIN.	-10
ABS RANGE MAX.	2.000
LINEARITY LIMIT	7.00
REAGENT BLANK	Yes
CONTAMINATING	No
DIFERENTIAL	No
FILTER 1	546
FILTER 2	None
TIMES MIX 1	0.0
INCUBATION 1	400
LAG PHASE	20
MIX 2	0.0
INCUBATION 2	0
MEASURE	80
MEASURE TYPE	Std
FACTOR	0.0
NORMAL RANGE	2.00 – 2.80

Vers. 0208