

RHEUMATOID FACTORS (RF)

Turbidimetry
LATEX

The application parameters comprised here constitute a guide to facilitate the validation of our reagents by the instrument. It is advisable to validate the use when there is any change in software or reagent versions.

Instrument: **AIRONE 200 (=AUTOHUMALYZER 900 S)**

Reagent preparation

Reagent 1: Use the Diluent.

Reagent 2: Use the Latex.

Instrument settings

UNITS	UI/mL
DECIMALS	1
REACTION TYPE	End Point
VOLUME SERUM	5
PLASMA	
URINE	
REAGENT 1	400
REAGENT 2	100
ABS RANGE MIN.	0.0
ABS RANGE MAX.	2.000
LINEARITY LIMIT	160
REAGENT BLANK	No
CONTAMINATING	Yes
DIFERENTIAL	No
FILTER 1	546
FILTER 2	None
TIMES MIX 1	1.0
INCUBATION 1	5
LAG PHASE	5
MIX 2	2.0
INCUBATION 2	150
MEASURE	1
MEASURE TYPE	Calibration curve
FACTOR	0.0
NORMAL RANGE	0.0 – 30.0

Vers. 0406