

# ALBUMIN

Spectrophotometric  
BROMOCRESOL GREEN

The application parameters comprised here constitute a guide to facilitate the validation of our reagents by the instrument. It is advisable to validate the use when there is any change in software or reagent versions.

## Instrument: **BS-300**

### Reagent preparation

Reagent is ready to use.

### Instrument settings

PARAMETERS							
N°	...	Test	<b>ALB</b>	Method	<b>End-Point</b>	Direction	<b>Ascend</b>
Unit	<b>g/L</b>	Decimal	<b>0</b>	Prim.Wavelength	<b>630</b>	Sec. Wavelength	<b>700</b>
Sample Volume	<b>3</b>	R1 Volume	<b>300</b>	R2 Volume	<b>0</b>	Line Limit	<b>N/A</b>
Incubation	--	Check Antigen	<b>No</b>	Response		Sample Volume	...
Reaction	<b>0 -- 10</b>	Substrate	<b>N/A</b>	Lower	--	Dilution	...
				Upper	--		
R1 Blank		Mix R Blank		Linearity		Full Name	...
Lower	...	Lower	...	Lower	<b>1.4</b>	Print Number	...
Upper	...	Upper	...	Upper	<b>70</b>		
Calibration							
Method	<b>Linear</b>	Std. 0	<b>0</b>				
		Std. 1	<b>(*)</b>				
		Std. 2	--				
		Std. 3	--				
		Std. 4	--				
		Std. 5	--				
						(*) Enter the standard value	
						... Value entered by the user	