

FERRITIN

Turbidimetry
LATEX

The application parameters comprised here constitute a guide to facilitate the validation of our reagents by the instrument. It is advisable to validate the use when there is any change in software or reagent versions.

Instrument: **BS-300**

Reagent preparation

Working Reagent: Pour the contents of a Reagent B vial (shake gently before using) into a Reagent A bottle. Mix thoroughly.
Stable for 20 days at 2-8°C.

Instrument settings

PARAMETERS							
N°	...	Test	FERR	Method	Fixed Time	Direction	Ascend
Unit	µg/L	Decimal	0	Prim.Wavelength	546	Sec. Wavelength	0
Sample Volume	6	R1 Volume	200	R2 Volume	0	Line Limit	20
Incubation	--	Check Antigen	No	Response		Sample Volume	...
Reaction	3 -- 28	Substrate	N/A	Lower	--	Dilution	...
				Upper	--		
R1 Blank		Mix R Blank		Linearity		Full Name	...
Lower	...	Lower	...	Lower	4	Print Number	...
Upper	...	Upper	...	Upper	500		
Calibration							
Method	Curve	Std. 0	0				
		Std. 1	(*) x 0,125				
		Std. 2	(*) x 0,25				
		Std. 3	(*) x 0,5				
		Std. 4	(*) x 0,75				
		Std. 5	(*) x 1				
						(*) Enter the standard value	
						... Value entered by the user	

Version 0707