

# ALKALINE PHOSPHATASE (ALP)

Continuous-spectrophotometric  
DEA BUFFER

The application parameters comprised here constitute a guide to facilitate the validation of our reagents by the instrument. It is advisable to validate the use when there is any change in software or reagent versions.

Instrument: **BT 3000 Plus (=ENVOY 500, FALCOR 350)**

## Reagent preparation

Working Reagent: Pour the contents of a Reagent B bottle into a Reagent A bottle. Mix gently.  
Stable for 2 months at 2-8°C.

## Instrument settings

Test Methodology		Units - Serum	U/L
Kind of Process	Kinetic	Units - Urine	
Type	Linear	Number of needle washes	2/2
Filters	405	Number of cuvette washes	2
Reaction Direction	Increase	Dynamic Blank	Inactive
Reagent #1 (µL)	300	Reagent Blank	Every Run
Reagent #2 (µL)	-	Reagent Limit (mABS)	1500
Sample Starter	Inactive	Curve Acceptance (%)	100
Delay Time (sec.)	0	Instrument Factor	0.000
Incubation Time (sec.)	60	Shift	0.000
Reading Time (sec.)	210		
- SERUM -		- URINE -	
Name	ALP DEA	Name	
Sample (µL)	6	Sample (µL)	
Pre-Dilution	1.00	Pre-Dilution	
<u>Dilution</u>		<u>Dilution</u>	
Factor	2.00	Factor	
Test Limit (Conc.)	690	Test Limit (Conc.)	
Initial ABS (mABS)	0	Initial ABS (mABS)	
Final ABS (mABS)	1900	Final ABS (mABS)	
Delta Max. ABS (mABS)	500	Delta Max. ABS (mABS)	
<u>Normal Values</u>		<u>Normal Values</u>	
Man	0 - 270	Man	
Woman	0 - 240	Woman	
Child		Child	
Calibrator			
Std. 1	0		
Std. 2	(*)		
		(*) Enter the assigned value of the calibrator	