

LIPASE

Color

The application parameters comprised here constitute a guide to facilitate the validation of our reagents by the instrument. It is advisable to validate the use when there is any change in software or reagent versions.

Instrument: **CLINLINE 150 (=LISA, MASCOTT Plus)**

Reagent preparation

Reagent 1: use the Reagent A.

Reagent 2: use the Reagent B.

Instrument settings

TEST NAME	LIPASE	BLANK = STD	
SHORT NAME	LIP	N° OF STD	
UNITS	U/L	STANDARD 1 VALUE	
ASSAY TYPE	KINETIC	STANDARD 1 POSITION	
FILTER VALUE	580	STANDARD 2 VALUE	
1 ST READING = 0	NO	STANDARD 2 POSITION	
LAG PHASE 1	12	N° REP. STANDARD/CONTROL	2
N° OF MEASURES	7	CONTROL VALUE	...
REAGENT 1 VOLUME	200	CONTROL POSITION	...
REAGENT 1 DILUTION	0	CONTROL DEVIATION	...
REAGENT 1 POSITION	...	PREDILUTION RATE	1
REAGENT 2 VOLUME	100	POSTDILUTION RATE	2
REAGENT 2 DILUTION	0	DILUENT	WATER
REAGENT 2 POSITION	...	RINSE TYPE	2
SAMPLE VOLUME	3	UP. NORMAL VALUE	48
SAMPLE DILUENT		LOW NORMAL VALUE	12
ACTIVATION	REAGENT 2	LINEARITY LIMIT	
LAG PHASE 2T	7	LOWER BLANK LIMIT	0
FACTOR		UPPER BLANK LIMIT	500
STD CALCULATION		LIMIT BLANK ACT.	5
		ODT 1 – ODT 0 LIMIT	500
... Data entered by the operator			
* Assigned value			