

α-1-ACID GLYCOPROTEIN

Turbidimetry

The application parameters comprised here constitute a guide to facilitate the validation of our reagents by the instrument. It is advisable to validate the use when there is any change in software or reagent versions.

Instrument: **CLINLINE 150 (=LISA, MASCOTT Plus)**

Reagent preparation

Reagent is ready to use.

Instrument settings

TEST NAME	α-1-Acid Glicoprotein	BLANK = STD	YES
SHORT NAME	AGP	N° OF STD	5
UNITS	mg/dL	STANDARD 1 VALUE	* x 0.125
ASSAY TYPE	E.P. STD	STANDARD 1 POSITION	...
FILTER VALUE	340	STANDARD 2 VALUE	* x 0.25
1 ST READING = 0	NO	STANDARD 2 POSITION	...
LAG PHASE 1	0	STANDARD 3 VALUE	* x 0.5
N° OF MEASURES	16	STANDARD 3 POSITION	...
REAGENT 1 VOLUME	400	STANDARD 4 VALUE	* x 0.75
REAGENT 1 DILUTION	0	STANDARD 4 POSITION	...
REAGENT 1 POSITION	...	STANDARD 5 VALUE	* x 1
REAGENT 2 VOLUME		STANDARD 5 POSITION	...
REAGENT 2 DILUTION		N° REP. STANDARD/CONTROL	1
REAGENT 2 POSITION		CONTROL VALUE	*
SAMPLE VOLUME	4	CONTROL POSITION	...
SAMPLE DILUENT		CONTROL DEVIATION	*
ACTIVATION	ANY	PREDILUTION RATE	1
LAG PHASE 2T		POSTDILUTION RATE	5
FACTOR		DILUENT	WATER
STD CALCULATION	4 DEG	RINSE TYPE	3
		UP. NORMAL VALUE	120
		LOW NORMAL VALUE	50
		LINEARITY LIMIT	170
		LOWER BLANK LIMIT	0
		UPPER BLANK LIMIT	1500

... Data entered by the operator

* Calibrator assigned value