

# RHEUMATOID FACTORS

Turbidimetry  
LATEX

The application parameters comprised here constitute a guide to facilitate the validation of our reagents by the instrument. It is advisable to validate the use when there is any change in software or reagent versions.

## Instrument: **CLINLINE 150 (=LISA, MASCOTT Plus)**

### Reagent preparation

Reagent 1: use the Diluent. Stable for 28 days at 2-8 °C.

Reagent 2: use the Latex. Shake the Latex vial gently before using.

### Instrument settings

TEST NAME		BLANK = STD	YES
	<b>RF</b>		
SHORT NAME	RF	N° OF STD	5
UNITS	U/mL	STANDARD 1 VALUE	* x 0.125
ASSAY TYPE	E.P. STD	STANDARD 1 POSITION	...
FILTER VALUE	620	STANDARD 2 VALUE	* x 0.25
1 <sup>ST</sup> READING = 0	NO	STANDARD 2 POSITION	...
LAG PHASE 1	0	STANDARD 3 VALUE	* x 0.5
N° OF MEASURES	4	STANDARD 3 POSITION	...
REAGENT 1 VOLUME	240	STANDARD 4 VALUE	* x 0.75
REAGENT 1 DILUTION	0	STANDARD 4 POSITION	...
REAGENT 1 POSITION	...	STANDARD 5 VALUE	* x 1
REAGENT 2 VOLUME	60	STANDARD 5 POSITION	...
REAGENT 2 DILUTION	0	N° REP. STANDARD/CONTROL	1
REAGENT 2 POSITION	...	CONTROL VALUE	*
SAMPLE VOLUME	3	CONTROL POSITION	...
SAMPLE DILUENT		CONTROL DEVIATION	*
ACTIVATION	ANY	PREDILUTION RATE	1
LAG PHASE 2T		POSTDILUTION RATE	5
FACTOR		DILUENT	WATER
STD CALCULATION	4 DEG	RINSE TYPE	3
		UP. NORMAL VALUE	20.0
		LOW NORMAL VALUE	0.0
		LINEARITY LIMIT	160
		LOWER BLANK LIMIT	0
		UPPER BLANK LIMIT	1500

... Data entered by the operator

\* Calibrator assigned value