

# ALANINE AMINOTRANSFERASE (ALT/GPT)

Continuous-spectrophotometric  
IFCC

The application parameters comprised here constitute a guide to facilitate the validation of our reagents by the instrument. It is advisable to validate the use when there is any change in software or reagent versions.

## Instruments: **CX-4/5/7**

### Reagent preparation

Working Reagent: Pour the contents of the Reagent B into the Reagent A bottle. Mix gently.  
Stable for 2 months at 2-8°C.

### Instrument settings

|                        |          |                            |                        |
|------------------------|----------|----------------------------|------------------------|
| Test Name              | ALT      | Calculation Factor         | 3333                   |
| Reaction Type          | RATE 1   | Math Model                 | LINEAR                 |
| Reaction Direction     | NEGATIVE | Cal. Time Limit            | 336 hr.                |
| Units                  | U/L      | N° of Calibrators          | 0                      |
| Decimal Precision      | X        | Secondary Wavelength       | 700                    |
| Primary Wavelength     | 340 nm   |                            |                        |
| Sample Volume          | 25 µL    | <b>CALIBRATION</b>         | <b>MULTIPOINT SPAN</b> |
| Primary Inject Reagent |          | N° 1                       |                        |
| A:                     | 250 µL   | N° 2                       |                        |
| B:                     |          | N° 3                       |                        |
| C:                     |          | N° 4                       |                        |
|                        |          | N° 5                       |                        |
|                        |          | N° 6                       |                        |
| <b>REAGENT BLANK</b>   |          | <b>REACTION</b>            |                        |
| Start Read             | 250      | Start Read                 | 120                    |
| End Read               | 300      | End Read                   | 180                    |
| Low ABS Limit          | -1.500   | Low ABS Limit              | -1.500                 |
| High ABS Limit         | 1.500    | High ABS Limit             | 1.500                  |
| <b>USABLE RANGE</b>    |          | <b>SUBSTRATE DEPLETION</b> |                        |
| Lower Limit            | 0        | Initial Rate               | -99.999                |
| Upper Limit            | 500      | Delta ABS                  | 1.500                  |
|                        |          | * assigned standard value  |                        |