

TRIGLYCERIDES

Enzymatic-spectrophotometric
GLYCEROL PHOSPHATE OXIDASE/PEROXIDASE

The application parameters comprised here constitute a guide to facilitate the validation of our reagents by the instrument. It is advisable to validate the use when there is any change in software or reagent versions.

Instruments: **CX-4/5/7**

Reagent preparation

Reagent is ready to be used.

Instrument settings

Test Name	TRIG	Calculation Factor	0
Reaction Type	ENDPOINT 2	Math Model	LINEAR
Reaction Direction	POSITIVE	Cal. Time Limit	336 hr.
Units	mg/dL	Nº of Calibrators	1
Decimal Precision	X	Secondary Wavelength	650
Primary Wavelength	520 nm		
Sample Volume	3 µL	CALIBRATION	MULTIPOINT SPAN
Primary Inject Reagent		Nº 1	*
A:	300 µL	Nº 2	
B:		Nº 3	
C:		Nº 4	
		Nº 5	
		Nº 6	
REAGENT BLANK		REACTION	
Start Read	250	Start Read	60
End Read	300	End Read	300
Low ABS Limit	-0.100	Low ABS Limit	0
High ABS Limit	0.300	High ABS Limit	2.000
USABLE RANGE		SUBSTRATE DEPLETION	
Lower Limit	1	Initial Rate	99.999
Upper Limit	600	Delta ABS	2.000
		* assigned standard value	