

# RHEUMATOID FACTORS

Turbidimetry  
LATEX

## Linear Calibration

The application parameters comprised here constitute a guide to facilitate the validation of our reagents by the instrument. It is advisable to validate the use when there is any change in software or reagent versions.

## Instruments: **COBAS MIRA / MIRA Plus**

### Reagent preparation

Reagent 1: use the Diluent. Stable for 28 days at 2-8 °C.

Reagent 2: use the Latex. Shake the Latex vial gently before using.

### Instrument settings

<b>GENERAL</b>		<b>CALCULATION</b>	
MEASUREMENT MODE	ABSORB	SAMPLE LIMIT	NO SPACE
REACTION MODE	R- S-SR1 3	REACTION DIRECTION	INCREASE 1
CALIBRATION MODE	CALIBRATOR 2	CHECK	ON 1
REAGENT BLANK	REAG/DIL 2	CONVERSION FACTOR	1
CLEANER	NO 1	OFFSET	0
WAVELENGTH	600 nm 5	TEST RANGE LOW	0
DECIMAL POSITION	1 (IU/mL)	HIGH	120 IU/mL
UNIT	27 (IU/mL)	NORM RANGE LOW	0 IU/mL
		HIGH	30 IU/mL
<b>ANALYSIS</b>		NUMBER OF STEPS	1
POST DIL. FACTOR	2	CALCULATION STEP A	ENDPOINT 1
CONC. FACTOR	NO SPACE	READING FIRST	CB
SAMPLE CYCLE	1	LAST	7
VOLUME	3 µL		
DILUTION NAME	H2O 0	<b>CALIBRATION</b>	
VOLUME	10 µL	CALIB. INTERVAL	ON REQUEST 3
REAGENT CYCLE	1	TIME	NO SPACE
VOLUME	240 µL	REAGENT RANGE LOW	NO SPACE
START R1 CYCLE	2	HIGH	NO SPACE
VOLUME	60 µL	BLANK RANGE LOW	NO SPACE
DILUTION NAME	H2O 0	HIGH	NO SPACE
VOLUME	10 µL		
		<b>STANDARD POS.</b>	(*)
		STD 1	(**)
(*) Enter the position of the CAL rack			
(**) Enter the value of the calibrator			