

LIPASE

Color

The application parameters comprised here constitute a guide to facilitate the validation of our reagents by the instrument. It is advisable to validate the use when there is any change in software or reagent versions.

Instrument: HITACHI 704

Reagent preparation

Reagent 1: use the Reagent A.

Reagent 2: use the Reagent B.

Instrument settings

TEST	LIPASE
ASSAY CODE	RATE A - 24 - 32
SAMPLE VOLUME (μL)	5
R1 VOLUME (μL)	330 - 20 - 0
R2 VOLUME (μL)	165 - 20 - 0
WAVELENGTH (nm)	700 - 570
CALIB. METHOD	1 - 0 - 0
STD 1 CONC. POS.	0 - 1
STD 2 CONC. POS.	(*) - 2
STD 3 CONC. POS.	0 - 0
STD 4 CONC. POS.	0 - 0
STD 5 CONC. POS.	0 - 0
STD 6 CONC. POS.	0 - 0
UNIT	U/L
SD LIMIT	0.1
DUPLICATE LIMIT	100
SENSITIVITY LIMIT	0
ABS. LIMIT	10000 - 0
PROZONE LIMIT	0 - 0
EXPECTED VALUE	12 - 48
INSTRUMENT FACTOR	1.0
(*) Input the assigned value of the calibrator	

Vers. 0605