

# TOTAL BILIRUBIN

Spectrophotometric  
DIAZOTIZED SULFANILIC

The application parameters comprised here constitute a guide to facilitate the validation of our reagents by the instrument. It is advisable to validate the use when there is any change in software or reagent versions.

## Instrument: **HITACHI 705**

### Reagent preparation

Working Reagent: Transfer the contents of one Reagent B-T vial into a Reagent A-T bottle. Mix thoroughly.  
Stable for 20 days at 2-8 °C.

### Instrument settings

TEST CODE	# BIL-T
ASSAY CODE	END POINT 1
SAMPLE VOLUME	20
R1 VOL. (µL)	320 - ...
R2 VOL. (µL)	80 - 1 - NO
WAVELENGTH 1	700
WAVELENGTH 2	546
RGT. BLK. ABS.	0
RGT. BLK. CONC.	0
STD. CONC.	*
FACTOR	...
STD. ABS. ALLOWANCE	10 %
NORMAL RANGE L	0.00 mg/dL
NORMAL RANGE H	1.10 mg/dL
ABS. LIMIT (RATE)	0
TEST CODE	# - #
	# Data entered by the operator
	* Assigned value

Vers. 0407