

CHOLINESTERASE

Continuous-spectrophotometric
BUTYRYLTHIOCHOLINE

The application parameters comprised here constitute a guide to facilitate the validation of our reagents by the instrument. It is advisable to validate the use when there is any change in software or reagent versions.

Instrument: HITACHI 705

Reagent preparation

Working Reagent: Pour the contents of the Reagent B into the Reagent A bottle. Mix gently.
Stable for 3 days at 2-8°C.

Instrument settings

TEST CODE	# (CHE)
ASSAY CODE	RATE 2 - 5 - 10
SAMPLE VOLUME	6
R1 VOL. (µL)	360 - ...
R2 VOL. (µL)	1 - 1 - NO
WAVELENGTH 1	505
WAVELENGTH 2	415
RGT. BLK. ABS.	0
RGT. BLK. CONC.	0
STD. CONC.	*
FACTOR	...
STD. ABS. ALLOWANCE	10 %
NORMAL RANGE L	3930 U/L
NORMAL RANGE H	11500 U/L
ABS. LIMIT (RATE)	4000
TEST CODE	#
	# Data entered by the operator
	* Assigned value

Vers. 0601