

# CREATINE KINASE

IFCC

The application parameters comprised here constitute a guide to facilitate the validation of our reagents by the instrument. It is advisable to validate the use when there is any change in software or reagent versions.

## Instrument: **HITACHI 705**

### Reagent preparation

Working Reagent: Pour the contents of the Reagent B into the Reagent A bottle. Mix gently.  
 Stable for 15 days at 2-8°C. Working Reagent must be protected from light.

### Instrument settings

TEST CODE	# (CK)
ASSAY CODE	RATE 2 - 24 - 28
SAMPLE VOLUME	16
R1 VOL. (µL)	350 - NO
R2 VOL. (µL)	1 - ... - NO
WAVELENGTH 1	376
WAVELENGTH 2	340
RGT. BLK. ABS.	0
RGT. BLK. CONC.	0
STD. CONC.	*
FACTOR	...
STD. ABS. ALLOWANCE	10 %
NORMAL RANGE L	0 U/L
NORMAL RANGE H	25 U/L
ABS. LIMIT (RATE)	11000
TEST CODE	#
	# Data entered by the operator
	* Assigned value

Vers. 0610