

CREATININE

Kinetic-spectrophotometric
ALKALINE PICRATE

The application parameters comprised here constitute a guide to facilitate the validation of our reagents by the instrument. It is advisable to validate the use when there is any change in software or reagent versions.

Instrument: HITACHI 705

Reagent preparation

Working Reagent: Mix equal volumes of Reagent A and Reagent B. Mix thoroughly.
Stable for 2 months at 2-8°C.

Instrument settings

TEST CODE	# (CREA)
ASSAY CODE	RATE 2 - 3 - 6
SAMPLE VOLUME	20
R1 VOL. (µL)	400 - ...
R2 VOL. (µL)	1 - 1 - NO
WAVELENGTH 1	570
WAVELENGTH 2	505
RGT. BLK. ABS.	0
RGT. BLK. CONC.	0
STD. CONC.	* - * - #
FACTOR	
STD. ABS. ALLOWANCE	10 %
NORMAL RANGE L	0.50 mg/dL
NORMAL RANGE H	1.10 mg/dL
ABS. LIMIT (RATE)	15000
TEST CODE	#
	# Data entered by the operator
	* Assigned value

Vers. 0312