

COMPLEMENT COMPONENT C4

Turbidimetry

The application parameters comprised here constitute a guide to facilitate the validation of our reagents by the instrument. It is advisable to validate the use when there is any change in software or reagent versions.

Instrument: **HITACHI 705**

Reagent preparation

Reagent is ready to use.

Instrument settings

| | |
|---------------------|--------------------------------|
| TEST CODE | # (C4) |
| ASSAY CODE | END POINT 1 |
| SAMPLE VOLUME | 10 |
| R1 VOL. (µL) | 400 - ... |
| R2 VOL. (µL) | 1 - 1 - NO |
| WAVELENGTH 1 | 700 |
| WAVELENGTH 2 | 340 |
| RGT. BLK. ABS. | 0 |
| RGT. BLK. CONC. | 0 |
| STD. CONC. | * - * - * |
| FACTOR | ... |
| STD. ABS. ALLOWANCE | 10 % |
| NORMAL RANGE L | 10 mg/dL |
| NORMAL RANGE H | 40 mg/dL |
| ABS. LIMIT (RATE) | 0 |
| TEST CODE | # |
| | # Data entered by the operator |
| | * Assigned value |

Vers. 0605