

FERRITIN

Turbidimetry
LATEX

The application parameters comprised here constitute a guide to facilitate the validation of our reagents by the instrument. It is advisable to validate the use when there is any change in software or reagent versions.

Instrument: **HITACHI 705**

Reagent preparation

Working Reagent: Mix four volumes of Diluent with one volume of Latex.

Stable for 1 day at 2-8°C.

Shake the latex vial gently before using.

Instrument settings

TEST CODE	# (FERRI)
ASSAY CODE	END POINT 1
SAMPLE VOLUME	20
R1 VOL. (µL)	400 - ...
R2 VOL. (µL)	- ...
WAVELENGTH 1	700
WAVELENGTH 2	546
RGT. BLK. ABS.	0
RGT. BLK. CONC.	0
STD. CONC.	* - * - *
FACTOR	...
STD. ABS. ALLOWANCE	10 %
NORMAL RANGE L	20 µg/L
NORMAL RANGE H	250 µg/L
ABS. LIMIT (RATE)	0
TEST CODE	#
	# Data entered by the operator
	* Assigned value

Vers. 0605