

# IMMUNOGLOBULIN G (IgG)

Turbidimetry

The application parameters comprised here constitute a guide to facilitate the validation of our reagents by the instrument. It is advisable to validate the use when there is any change in software or reagent versions.

## Instrument: **HITACHI 705**

### Reagent preparation

Reagent is ready to use.

### Instrument settings

|                     |                                |
|---------------------|--------------------------------|
| TEST CODE           | # (IG G)                       |
| ASSAY CODE          | END POINT 1                    |
| SAMPLE VOLUME       | 5                              |
| R1 VOL. (µL)        | 500 - ...                      |
| R2 VOL. (µL)        | 1 - 1 - NO                     |
| WAVELENGTH 1        | 700                            |
| WAVELENGTH 2        | 546                            |
| RGT. BLK. ABS.      | 0                              |
| RGT. BLK. CONC.     | 0                              |
| STD. CONC.          | * - * - *                      |
| FACTOR              | ...                            |
| STD. ABS. ALLOWANCE | 10 %                           |
| NORMAL RANGE L      | 700 mg/dL                      |
| NORMAL RANGE H      | 1600 mg/dL                     |
| ABS. LIMIT (RATE)   | 0                              |
| TEST CODE           | #                              |
|                     | # Data entered by the operator |
|                     | * Assigned value               |

Vers. 0605