

# CALCIUM

Spectrophotometric  
ARSENAZO III

The application parameters comprised here constitute a guide to facilitate the validation of our reagents by the instrument. It is advisable to validate the use when there is any change in software or reagent versions.

## Instrument: **HITACHI 902 (=HITACHI 7020)**

### Reagent preparation

Reagent is ready to be used.

### Instrument settings

|                                  |         |                          |             |
|----------------------------------|---------|--------------------------|-------------|
| 1 Test Name                      | CA      | 29 Calib. Conc. 4        | 0           |
| 2 Assay Code (Mthd)              | 1 Point | 30 Calib. Position 4     | 0           |
| 3 Assay Code (2. Test)           | 0       | 31 Calib. Conc. 5        | 0           |
| 4 Reaction Time                  | 3       | 32 Calib. Position 5     | 0           |
| 5 Assay Point 1                  | 8       | 33 Calib. Conc. 6        | 0           |
| 6 Assay Point 2                  | 0       | 34 Calib. Position 6     | 0           |
| 7 Assay Point 3                  | 0       | 35 S1 ABS                |             |
| 8 Assay Point 4                  | 0       | 36 K Factor              |             |
| 9 Wavelength (SUB)               |         | 37 K2 Factor             | 0           |
| 10 Wavelength (MAIN)             | 660     | 38 K3 Factor             | 0           |
| 11 Sample Volume                 | 4       | 39 K4 Factor             | 0           |
| 12 R1 Volume                     | 270     | 40 K5 Factor             | 0           |
| 13 R1 Position                   | ...     | 41 A Factor              | 0           |
| 14 R1 Bottle size                | Small   | 42 B Factor              | 0           |
| 15 R2 Volume                     | 0       | 43 C Factor              | 0           |
| 16 R2 Position                   | 0       | 44 SD Limit              | 0.1         |
| 17 R2 Bottle size                | Small   | 45 Duplicate Limit       | 200         |
| 18 R3 Volume                     | 0       | 46 Sensitivity Limit     | 0           |
| 19 R3 Position                   | 0       | 47 S1 ABS. Limit (L)     | -32000      |
| 20 R3 Bottle size                | Small   | 48 S1 ABS. Limit (H)     | 32000       |
| 21 Calib. Type (Type)            | Linear  | 49 ABS. Limit            | 0           |
| 22 Calib. Type (Weight)          | 0       | 50 ABS. Limit (D/I)      | Increase    |
| 23 Calib. Conc. 1                | 0.0     | 51 Prozone Limit         | 0           |
| 24 Calib. Position 1             | ...     | 52 Prozone Limit (U/D)   | Lower limit |
| 25 Calib. Conc. 2                | *       | 53 Prozone (End Point)   | 35          |
| 26 Calib. Position 2             | ...     | 54 Expected Value (L)    | 9.0         |
| 27 Calib. Conc. 3                | 0       | 55 Expected Value (H)    | 10.7        |
| 28 Calib. Position 3             | 0       | 56 Instrument Factor (a) | 1.0         |
|                                  |         | 57 Instrument Factor (b) | 0.0         |
| ... Data entered by the operator |         | 58 Key Setting           | ...         |
| * Assigned value                 |         |                          |             |