

GLUCOSE

Enzymatic-spectrophotometric
GLUCOSE OXIDASE/PEROXIDASE

The application parameters comprised here constitute a guide to facilitate the validation of our reagents by the instrument. It is advisable to validate the use when there is any change in software or reagent versions.

Instrument: **ILAB 300 / 300 Plus (=LIASYS)**

Reagent preparation

Reagent is ready to be used.

Instrument settings

| | | | | | | | | | | | |
|-------------------------------------|-------|-------------------------------|-----------------|----------|--------|---------------|--------|------|-----------|---------------------|----------------------|
| Description | | GLU | Reference Range | | | | | | | | |
| Unit | mg/dL | Male | NORMAL VALUES | | | HIGH VALUES | | | | | |
| Decimals | 0 | Female | 0 | 0 | 0 | 70 | 105 | ... | ... | ... | |
| LIS Code | ... | Children | 0 | 0 | 0 | 70 | 105 | ... | ... | ... | |
| Unit Factor | 1.00 | Alert | Low Alert | Very Low | Low | Normal Values | | High | Very High | High Alert | |
| Slope | 1.00 | Rerun | No | No | | | | | No | No | |
| Intercept | 0.00 | | | | | | | | | | |
| Reaction Type | | End Point | Parameters | | | | | | | | |
| Direction | Up | Times(sec) | Predilut. | S+R 1 | Reag 2 | Reag 3 | Incub. | Read | | | |
| E.P.Limit (Abs) | 1.0 | Dil./Rgt Code | | 0 | 0 | 0 | 250 | | | | |
| Depl.limit (Abs) | N/A | Lot Number | | GLU | 0 | | | | | | |
| First Limit (Abs) | N/A | Ratio/Vol (uL) | 1/1 | 300 | 0 | 0 | | | | | |
| Linear Factor | N/A | Rinse (uL) | | 0 | 0 | 0 | | | | | |
| Fit: | N/A | Sample (uL) | | 3 | | | | | | | |
| Lin Limit Low | 0 | RBL min (abs) | 0.0 | | | | | | | | |
| Lin Limit High | 500 | Max (abs) | 1.0 | | | | | | | | |
| Rerun when over | No | | | | | | | | | Filter 1 (nm) : 510 | |
| Calculation Model | | Standard | | | | | | | | | Filter 2 (nm) : none |
| Factor | N/A | RBL stability (days) : | 1 | | | | | | | | Bicr.Factor : 1 |
| Sample Blank | No | Calibration Stability (days): | 99 | | | | | | | | |
| | | Dynamic Controls (min): | none | | | | | | | | |
| (...) Value entered by the operator | | | | | | | | | | | |