

TOTAL BILIRUBIN

Spectrophotometric
DIAZOTIZED SULFANILIC

The application parameters comprised here constitute a guide to facilitate the validation of our reagents by the instrument. It is advisable to validate the use when there is any change in software or reagent versions.

Instrument: **ILAB 900 (=SHIMADZU CL-7200/7300)**

Reagent preparation

Working Reagent: Transfer the contents of one Reagent B-T vial into a Reagent A-T bottle. Mix thoroughly.
Stable for 20 days at 2-8 °C.

Instrument settings

TEST N°:		CURVE N°:	
TEST NAME	BIL-T	NORMAL RANGE MALE L - H	0.0 – 1.1
BAR CODE		FEMALE L - H	0.0 – 1.1
MEASURING METHOD	Endpoint	VALID RANGE LL - HH	0 – 15
CALIB. METHOD	1 point	CORRECTION CONST. 1,2	1 – 0
UNIT	mg/dL	W1 – ABS INC / DEC, LIMIT	Inc – 2000
MEASURING POINT	43	HEMOLYSIS LIMIT	-
PHOTOMETRIC METHOD	Double	ICTERIC LIMIT	-
WAVELENGTHS 1, 2, 3	540 – 700	LIPEMIC LIMIT	-
SAMPLE VOL. N, D, I	30 – 15 – 30	CALIB. REPLICATES	5
R1 V. – STIR MODE (Nor)	280 – Yes – 0	REPLICATE, SENS. LIMITS	-
R2 V. – STIR MODE		SLOPE – OFFSET LIMITS	-
N°. DECIMAL POINTS	2	PARTNER TEST, CONST.	-
SAMPLE TYPE	Serum	PROZONE INC / DEC, LIMIT	-
... Data entered by the operator		* assigned value	