

# APOLIPOPROTEIN B (Apo B)

Turbidimetry

The application parameters comprised here constitute a guide to facilitate the validation of our reagents by the instrument. It is advisable to validate the use when there is any change in software or reagent versions.

Instrument: **ILAB 900 (=SHIMADZU CL-7200/7300)**

### Reagent preparation

Reagent 1: Use the Reagent A.

Reagent 2: Use the Reagent B.

### Instrument settings

<b>TEST Nº.:</b> TEST NAME  <div style="text-align: center;"><b>APO B</b></div> BAR CODE MEASURING METHOD CALIB. METHOD UNIT MEASURING POINT PHOTOMETRIC METHOD WAVELENGTHS 1, 2, 3 SAMPLE VOL. N, D, I R1 V. – STIR MODE (Nor) R2 V. – STIR MODE (Nor) Nº. DECIMAL POINTS SAMPLE TYPE  ... Data entered by the operator	<b>CURVE Nº.:</b> NORMAL RANGE MALE L - H  <div style="text-align: center;">FEMALE L - H</div> VALID RANGE LL - HH CORRECTION CONST. 1,2 W1 – ABS INC / DEC, LIMIT HEMOLYSIS LIMIT ICTERIC LIMIT LIPEMIC LIMIT CALIB. REPLICATES REPLICATE, SENS. LIMITS SLOPE – OFFSET LIMITS PARTNER TEST, CONST. PROZONE INC / DEC, LIMIT  * assigned value
	63 – 133  63 – 133 2 – 300 1 – 0 Inc – 3000 - - - 5 2500 – 2500 999 – 2500 - -