

ANTI-STREPTOLYSIN O (ASO)

Turbidimetry
LATEX

The application parameters comprised here constitute a guide to facilitate the validation of our reagents by the instrument. It is advisable to validate the use when there is any change in software or reagent versions.

Instrument: **ILAB 900 (=SHIMADZU CL-7200/7300)**

Reagent preparation

Reagent 1: use the Diluent. Stable for 20 days at 2-8 °C.

Reagent 2: use the Latex. Shake the Latex before pipetting. Stable for 20 days at 2-8 °C.

Instrument settings

TEST Nº:		CURVE Nº:	
TEST NAME		NORMAL RANGE MALE L - H	0 - 200
	ASO		
BAR CODE		FEMALE L - H	0 - 200
MEASURING METHOD	Endpoint	VALID RANGE LL - HH	0 - 800
CALIB. METHOD	1 point	CORRECTION CONST. 1,2	1 - 0
UNIT	IU/mL	W1 - ABS INC / DEC, LIMIT	Inc - 2500
MEASURING POINT	16 - 26	HEMOLYSIS LIMIT	-
PHOTOMETRIC METHOD	Single	ICTERIC LIMIT	-
WAVELENGTHS 1, 2, 3	540	LIPEMIC LIMIT	-
SAMPLE VOL. N, D, I	3 - 2 - 30	CALIB. REPLICATES	5
R1 V. - STIR MODE (Dil)	240 - Yes - 10	REPLICATE, SENS. LIMITS	2500 - 2500
R2 V. - STIR MODE (Dil)	60 - Yes - 10	SLOPE - OFFSET LIMITS	999 - 2500
Nº. DECIMAL POINTS	0	PARTNER TEST, CONST.	-
SAMPLE TYPE	Serum	PROZONE INC / DEC, LIMIT	-
... Data entered by the operator		* assigned value	

Version 0410