

ALBUMIN

Spectrophotometric
BROMOCRESOL GREEN

The application parameters comprised here constitute a guide to facilitate the validation of our reagents by the instrument. It is advisable to validate the use when there is any change in software or reagent versions.

Instrument: RA-XT (=OPERA)

Reagent preparation

Reagent is ready to use.

Instrument settings

CHEM. #	(*)
NAME	ALB
IMMUNOASSAY	NO
TYPE	ENDPOINT
INVERSE	NO
% SAMPLE VOL.	6 (3 µL)
FILTER P.	600
BIC. CHEM.	NO
DELAY TIME	1:00
BLANK	NO
% REAGENT VOL.	60 (300 µL)
2 ND REAGENT	NO
UNITS	g/L
UNIT FACTOR	1.000
DECIMAL	1
RBL. LOW	0.000
RBL. HIGH	0.450
RANGE LOW	0.0
RANGE HIGH	70.0
CAL. FACTOR	(**)
STANDARD VALUE	(***)
REAGENT RATE	0.000
NORMAL LOW	30
NORMAL HIGH	50
SLOPE	1.00
INTERCEPTION	0.00
ENDPOINT LIMIT	0.030
	(*) Data entered by the operator
	(**) CAL. FACTOR is determined by a calibration assay
	(***) Assigned value of the standard