

CHOLINESTERASE

Continuous-spectrophotometric
BUTYRYLTHIOCHOLINE

The application parameters comprised here constitute a guide to facilitate the validation of our reagents by the instrument. It is advisable to validate the use when there is any change in software or reagent versions.

Instrument: RA-XT (=OPERA)

Reagent preparation

Working Reagent: Pour the contents of the Reagent B into the Reagent A bottle. Mix gently.
Stable for 3 days at 2-8°C.

Instrument settings

CHEM. #	(*)
NAME	CHE
IMMUNOASSAY	NO
TYPE	0
INVERSE	YES
% SAMPLE VOL.	10
FILTER P.	405
BIC. CHEM.	NO
DELAY TIME	2:00
BLANK	NO
% REAGENT VOL.	60
2 ND REAGENT	NO
UNITS	U/L
UNIT FACTOR	1.000
DECIMAL	0
RBL. LOW	0.0
RBL. HIGH	0.600
RANGE LOW	0
RANGE HIGH	25000
CAL. FACTOR (**)	65804
REAGENT RATE	0.000
NORMAL LOW	3930
NORMAL HIGH	11500
SLOPE	1.00
INTERCEPTION	0
C1*10E-6	0.00
C2*10E-6	99999
D1*10E-6	10.00
DELTA#	0.010
1=ACCEPT	1 (YES)
	(*) Data entered by the operator (**) CAL. FACTOR is determined by a calibration assay