

FRUCTOSAMINE

Kinetic-spectrophotometric
NTB

The application parameters comprised here constitute a guide to facilitate the validation of our reagents by the instrument. It is advisable to validate the use when there is any change in software or reagent versions.

Instrument: RA-XT (=OPERA)

Reagent preparation

Reagent is ready to be used.

Instrument settings

CHEM. #	(*)
NAME	FRUC
IMMUNOASSAY	NO
TYPE	FIRST ORDER
INVERSE	NO
% SAMPLE VOL.	30 (15 µL)
FILTER P.	550
BIC. CHEM.	NO
DELAY TIME	7:00
INCUBATE	2:30
BLANK TYPE	NO
% REAGENT VOL.	60 (300 µL)
2 ND REAGENT	NO
UNITS	mmol/L
UNIT FACTOR	1.000
DECIMAL	2
RBL. LOW	0.000
RBL. HIGH	0.300
RANGE LOW	0.00
RANGE HIGH	7.00
CAL. FACTOR	(**)
STANDARD VALUE	(***)
REAGENT RATE	0.000
NORMAL LOW	2.00
NORMAL HIGH	2.80
SLOPE	1.00
INTERCEPTION	0.00
LINEAR FACTOR	2.65
ENDPOINT LIMIT	0.100
	(*) Data entered by the operator
	(**) CAL. FACTOR is determined by a calibration assay
	(***) Assigned value of the standard

Version 0301