

GAMMA-GLUTAMYLTRANSFERASE

Continuous-spectrophotometric
IFCC

The application parameters comprised here constitute a guide to facilitate the validation of our reagents by the instrument. It is advisable to validate the use when there is any change in software or reagent versions.

Instrument: **RA-XT (=OPERA)**

Reagent preparation

Working Reagent: Pour the contents of the Reagent B into the Reagent A bottle. Mix gently.
Stable for 2 months at 2-8 °C.

Instrument settings

CHEM. #	(*)
NAME	GGT
IMMUNOASSAY	NO
TYPE	ZERO ORDER
INVERSE	NO
% SAMPLE VOL.	60 (30 µL)
FILTER P.	405
BIC. CHEM.	NO
DELAY TIME	0:30
BLANK	NO
% REAGENT VOL.	60 (300 µL)
2 ND REAGENT	NO
UNITS	U/L
UNIT FACTOR	1.000
DECIMAL	0
RBL. LOW	0.150
RBL. HIGH	0.800
RANGE LOW	0.0
RANGE HIGH	300
CAL. FACTOR (**)	1587
REAGENT RATE	0.000
NORMAL LOW	15
NORMAL HIGH	86
SLOPE	1.00
INTERCEPTION	0:00
C1*10E-6	0.00
C2*10E-6	99999
D1*10E-6	5.60
DELTA#	0.010
	(*) Data entered by the operator
	(**) CAL. FACTOR is determined by a calibration assay

Version 0205