

# IRON

Spectrophotometric  
CHROMAZUROL B

The application parameters comprised here constitute a guide to facilitate the validation of our reagents by the instrument. It is advisable to validate the use when there is any change in software or reagent versions.

## Instrument: RA-XT (=OPERA)

### Reagent preparation

Reagent is ready to be used.

### Instrument settings

CHEM. #	(*)
NAME	IRON C
IMMUNOASSAY	NO
TYPE	ENDPOINT
INVERSE	NO
% SAMPLE VOL.	30 (15 µL)
FILTER P.	600
BIC. CHEM.	NO
DELAY TIME	4:00
BLANK TYPE	NO
% REAGENT VOL.	60 (300 µL)
2 <sup>ND</sup> REAGENT	NO
UNITS	µg/dL
UNIT FACTOR	1.000
DECIMAL	0
RBL. LOW	0.000
RBL. HIGH	0.050
RANGE LOW	0.00
RANGE HIGH	500
CAL. FACTOR	(**)
STANDARD VALUE	(***)
NORMAL LOW	55
NORMAL HIGH	155
SLOPE	1.00
INTERCEPTION	0.00
ENDPOINT LIMIT	0.020
	(*) Data entered by the operator
	(**) CAL. FACTOR is determined by a calibration assay
	(***) Assigned value of the standard

Version 0208