

LIPASE

Color

The application parameters comprised here constitute a guide to facilitate the validation of our reagents by the instrument. It is advisable to validate the use when there is any change in software or reagent versions.

Instrument: **RA-XT (=OPERA)**

Reagent preparation

Reagent 1: use the Reagent A.
 Reagent 2: use the Reagent B.

Instrument settings

CHEM. #	(*)
NAME	LIP
IMMUNOASSAY	NO
TYPE	ZERO ORDER
INVERSE	NO
% SAMPLE VOL.	6 (3 µL)
FILTER P.	550
BIC. CHEM.	NO
DELAY TIME	2:45
% REAGENT VOL.	40 (200 µL)
2 ND REAGENT	YES
% 2 ND REAGENT	20 (100 µL)
A2 DELAY	2:45
UNITS	U/L
UNIT FACTOR	1.000
DECIMAL	0
RBL. LOW	0.000
RBL. HIGH	0.800
RANGE LOW	3.0
RANGE HIGH	400
CAL. FACTOR (**)	
REAGENT RATE	0.000
NORMAL LOW	12
NORMAL HIGH	48
SLOPE	1.00
INTERCEPTION	0:00
C1*10E-6	0.00
C2*10E-6	99999
D1*10E-6	6.00
DELTA#	0.010
	(*) Data entered by the operator (**) CAL. FACTOR is determined by a calibration assay