

PROTEIN (URINE)

Spectrophotometric
PYROGALLOL RED

The application parameters comprised here constitute a guide to facilitate the validation of our reagents by the instrument. It is advisable to validate the use when there is any change in software or reagent versions.

Instrument: **RA-XT (=OPERA)**

Reagent preparation

Reagent is ready to be used.

Instrument settings

CHEM. #	(*)
NAME	PROT-UR
IMMUNOASSAY	NO
TYPE	ENDPOINT
INVERSE	NO
% SAMPLE VOL.	24 (12 µL)
FILTER P.	600
BIC. CHEMISTRY	YES
BIC. WAVELENGTH	550
BIC. TYPE	FIXED
BIC. FACTOR	1
DELAY TIME	3:00
BLANK	NO
% REAGENT VOL.	60 (300 µL)
2 ND REAGENT	NO
UNITS	mg/dL
UNIT FACTOR	1.000
DECIMAL	0
RBL. LOW	0.000
RBL. HIGH	0.150
RANGE LOW	0.0
RANGE HIGH	200
CAL. FACTOR	(**)
STANDARD VALUE	(***)
NORMAL LOW	3.0
NORMAL HIGH	14.0
SLOPE	1.00
INTERCEPTION	0.00
ENDPOINT LIMIT	0.50
	(*) Data entered by the operator
	(**) CAL. FACTOR is determined by a calibration assay
	(***) Assigned value of the standard