

TRIGLYCERIDES

Enzymatic-spectrophotometric
GLYCEROL PHOSPHATE OXIDASE/PEROXIDASE

The application parameters comprised here constitute a guide to facilitate the validation of our reagents by the instrument. It is advisable to validate the use when there is any change in software or reagent versions.

Instrument: **RA-XT (=OPERA)**

Reagent preparation

Reagent is ready to be used.

Instrument settings

| | |
|-------------------------|---|
| CHEM. # | (*) |
| NAME | TRIG |
| IMMUNOASSAY | NO |
| TYPE | ENDPOINT |
| INVERSE | NO |
| % SAMPLE VOL. | 6 (3 μ L) |
| FILTER P. | 500 |
| BIC. CHEM. | NO |
| DELAY TIME | 4:00 |
| BLANK TYPE | NO |
| % REAGENT VOL. | 60 (300 μ L) |
| 2 ND REAGENT | NO |
| UNITS | mg/dL |
| UNIT FACTOR | 1.000 |
| DECIMAL | 0 |
| RBL. LOW | 0.000 |
| RBL. HIGH | 0.050 |
| RANGE LOW | 0.00 |
| RANGE HIGH | 600 |
| CAL. FACTOR | (**) |
| STANDARD VALUE | (***) |
| NORMAL LOW | 60 |
| NORMAL HIGH | 150 |
| SLOPE | 1.00 |
| INTERCEPTION | 0.00 |
| ENDPOINT LIMIT | 0.020 |
| | (*) Data entered by the operator |
| | (**) CAL. FACTOR is determined by a calibration assay |
| | (***) Assigned value of the standard |