

C-REACTIVE PROTEIN (CRP)

Turbidimetry
LATEX

The application parameters comprised here constitute a guide to facilitate the validation of our reagents by the instrument. It is advisable to validate the use when there is any change in software or reagent versions.

Instrument: **RA-XT (=OPERA)**

Reagent preparation

Reagent 1: Use the Reagent A.

Reagent 2: Use the Reagent B.

Instrument settings

CHEM. #	(*)
NAME	CRP
IMMUNOASSAY	NO
TYPE	ENDPOINT
INVERSE	NO
% SAMPLE VOL.	5 (2.5 µL)
FILTER P.	550
BIC. CHEMISTRY	NO
DELAY TIME	2:00
% REAGENT VOL.	60 (300 µL)
2 ND REAGENT	YES
% 2 ND REAGENT	15 (75 µL)
A2 DELAY	2:00
UNITS	mg/dL
UNIT FACTOR	1.000
DECIMAL	0
RBL. LOW	0.000
RBL. HIGH	0.800
RANGE LOW	0.00
RANGE HIGH	150
CAL. FACTOR	(**)
STANDARD VALUE	(***)
NORMAL LOW	0.0
NORMAL HIGH	5.0
SLOPE	1.00
INTERCEPTION	0.00
ENDPOINT LIMIT	0.150
	(*) Data entered by the operator
	(**) CAL. FACTOR is determined by a calibration assay
	(***) Assigned value of the standard