

ALBUMIN

Spectrophotometric
BROMOCRESOL GREEN

The application parameters comprised here constitute a guide to facilitate the validation of our reagents by the instrument. It is advisable to validate the use when there is any change in software or reagent versions.

Instrument: **RX-DAYTONA (=BQA, CHEMIX-180)**

Reagent preparation

Reagent is ready to use.

Instrument settings

CHEMISTRY PARAMETERS				ALB	R1	Reagent Name	ALB	Volume 300 µL
Name				ALB	R1	Reagent Name	ALB	Volume 300 µL
Unit				g/L	R2	Reagent Name		
Assay Type				Endpoint				
Measuring Points	1	Start			Decimal points	0		
		End		5	Normal Range (Male)	35 - 50		
	2	Start			(Female)	35 - 50		
		End			(Child)	28 - 54		
Wavelength					Technical Range (Conc.)	1.4 - 70		
Prim.	660		Sec.		(mAbs/10)	-35000 - 35000		
Sampling Volume	3 µL				RPT Wash (R1)	Sys Water		Wash Sol.-3
Rerun	2 µL				(R2)	Sys Water		Wash Sol.-3
CALIBRATION				S1	0	Instrument Factor	a 1	b 0
Calculation	Linear			S2	...	Stirring Speed	R1 Mid	R2 Mid
... Data entered by the operator								

Vers. 0510