

URIC ACID

Enzymatic-spectrophotometric
URICASE/PEROXIDASE

The application parameters comprised here constitute a guide to facilitate the validation of our reagents by the instrument. It is advisable to validate the use when there is any change in software or reagent versions.

Instrument: SPECTRUM CCX

Reagent preparation

Reagent is ready to be used.

Instrument settings

TEST DEFINITION							
ENTRY NAME	URIC	SAMPLE (μL)	NORMAL	5.0			
REPORT NAME	URIC ACID		LOW				
RATIO REF.	URIC		HIGH	2.5			
TEST NUMBER	***	UNITS PRIM		mg/dL			
TEST TYPE	CALIBRATED	SEC. UNIT FACTOR					
MATH	1 PT CAL FACTOR	PRINT DIGITS		0			
REACTION DIRECTION	UP						
REAGENTS	1	INST MUL	1.000	INT	0.000		
TEMPERATURE	37	NORMAL (C)		2.4 TO 7.0			
TEST BLANK TYPE	REAGENT BLANK	SAMPLE DISP. DELAY (sec.)		0			
LINEAR MODEL CALIBRATION DEFINITION							
TEST NAME	URIC	TEST TYPE		CALIBRATED			
COMB. TEST		MATH MODEL		LIN REG	END PT		
CAL MODE	CAL ON CMD	CAL INTERVAL (hr.)		720			
INTCPT TOL (C)	-1000.00 TO 1000.00	REF CAL FACTOR		*			
% TOL OF CAL FACTOR	10	% TOL OF CAL		10			
CAL LEVEL	0						
CALIBRATOR	LEVEL (C)	REPLICATES					
WATER	0.000	1					
MCC1	***	1					
MCC2	***	1					
MCC3	***	1					
TEST PARAMETER FILE: REAGENT DEFINITION							
REAGENT NAME	URIC	REAGENT NUMBER	1 FOR TEST URIC ACID	LINEARITY [C]	1.000 TO 25.0		
REAGENT VOL (μL) NORM	236			INITIAL Ad	0.2		
LOW	0.00			ABS LIMIT (Ad)	0.7		
HIGH	0.00			REAGENT BLANK			
FIRST READ TIME (sec)	240			BEFORE WASH CYCLES	18		
LAST READ TIME (sec)	240			AFTER WASH CYCLES	18		
NUMBER OF READS	1			MIX TIME (sec)	1:00		
READ INTERVAL (sec)	60			COOLING	YES		
AUX REAG DISP (sec)	0			CONSTANT INTERCEPT	0.000		
RSM	0			SPECTRAL CORRECTION			
PRIM	SEC	USE IN	CONST	E.F.	LOW	HIGH	
500	/ 604	A	1.00	0.00	0.00	0.00	