

# ALBUMIN

Spectrophotometric  
BROMOCRESOL GREEN

The application parameters comprised here constitute a guide to facilitate the validation of our reagents by the instrument. It is advisable to validate the use when there is any change in software or reagent versions.

## Instruments: **BT TARGA 2000/3000 (=ATAC 8000)**

### Reagent preparation

Reagent is ready to use.

### Instrument settings

Test Type		<i>End Point Single</i>
Serum Starter		Inactive
Filters	(A/B)	630
Units		g/L
Test Method		With factor
Test Methodology		
Number of Washes		1
Delay Time	(Sec)	0
Inc. Reagent	(Sec)	60
Reading Time	(Sec)	20
Test Limit	(Conc)	70
Max Abs Delta	(mABS)	0.00
Reagent mABS Limit		800
Reagent A/B	( $\mu$ L)	300
Reaction Direction		Increasing
Reagent Dilution		1:1
Initial ABS	(mABS)	0.00
Curve Acceptance		20 %
Automatic Profile	(Inactive)	Inactive
Rerun Test Rgt BLK H:M		00:00
<b>Serum Parameters</b>		
Test's Name		Albumin
Sample Volume	( $\mu$ L)	3
Dilution Ratio		1:2
Min. Max. M		35 – 50
Min. Max. F		35 – 50
Min. Max. B		
(*) Input the assigned value of the calibrator		