

CREATININE

Kinetic-spectrophotometric
ALKALINE PICRATE

The application parameters comprised here constitute a guide to facilitate the validation of our reagents by the instrument. It is advisable to validate the use when there is any change in software or reagent versions.

Instrument: XL-600 (=XL-300+ISE)

Reagent preparation

Reagent 1: Use the Reagent A.

Reagent 2: Use the Reagent B.

Instrument settings

Test Code Test	CREA						Reported Name	CREATININE			
Assay Type	2 POINT							M1Start	M1End	M2Start	M2End
Wavelength	Primary	505		Secondary	700		Assay Points	19	19	27	27
							Con. Interval	*			
							Sample Repli.	1			
	Serum			Urine							
	Sample	Predil	Diluent	Sample	Predil	Diluent					
S. Vol. Normal	22	0	0	4	22	196		Vol.	Pos.	Size	
S. Vol. Decr	10	0	0	4	10	196	R 1	110	*	L	
S. Vol. Incr	30	0	0	4	30	196	R 2	110	*	L	
Std. Volume	22						Reagent Stability		Effective Days		
ABS Limit	0			0				Min		Max	
React. Dir.	<input type="radio"/>		Decr	<input checked="" type="radio"/>		Incr					
Prozone Limit	0		<input type="radio"/>	Upper	<input type="radio"/>	Lower	Reagent ABS	0		0	
Unit	mg/dL						Decimal Point	1			
							Tech. Serum Lim.	0		20	
							Tech. Urine Lim.	0		0	
							Panic Limit	*		*	
Normal Values	AGE	Male		Female							
		Min	Max	Max	Max						
Serum	Default	0.9	1.3	0.6	1.1		Auto Dil.	<input checked="" type="radio"/>	Yes	<input type="radio"/>	No
Serum							Y=aX+b	a =	1	b =	0
Urine values											
Calibration curve	Straight										

* Data entered by the operator

In the **Std. Volume (Pre/Norm/Dil)** field, enter the **S. Vol. Normal ((Pre/Norm/Dil)** values