

# PHOSPHORUS

Spectrophotometric  
PHOSPHOMOLYBDATE/UV

The application parameters comprised here constitute a guide to facilitate the validation of our reagents by the instrument. It is advisable to validate the use when there is any change in software or reagent versions.

## Instrument: **XL-600 (=XL-300+ISE)**

### Reagent preparation

Working Reagent: Mix 35 mL Reagent A + 15 mL Reagent B. Mix thoroughly. Stable for 12 months at 15-30 °C.

### Instrument settings

Test Code Test	<b>PHOS</b>		Reported Name				<b>PHOSPHORUS</b>					
Assay Type	<b>2 POINT</b>		M1Start	M1End	M2Start	M2End						
Wavelength	Primary	<b>340</b>	Secondary	<b>0</b>	Assay Points		<b>12</b>	<b>12</b>	<b>51</b>	<b>51</b>		
					Con. Interval		*					
					Sample Repli.		1					
	Serum			Urine								
	Sample	Predil	Diluent	Sample	Predil	Diluent	Vol.	Pos.	Size			
S. Vol. Normal	<b>3</b>	<b>0</b>	<b>0</b>	<b>20</b>	<b>3</b>	<b>180</b>	R 1	<b>220</b>	*	<b>L</b>		
S. Vol. Decr	<b>30</b>	<b>9</b>	<b>225</b>	<b>10</b>	<b>3</b>	<b>190</b>	R 2	<b>85</b>	*	<b>S</b>		
S. Vol. Incr	<b>6</b>	<b>0</b>	<b>0</b>	<b>40</b>	<b>3</b>	<b>160</b>	Reagent Stability		Effective Days			
Std. Volume	<b>3</b>											
ABS Limit	<b>0</b>				<b>0</b>							
React. Dir.	<input type="radio"/>		Decr	<input checked="" type="radio"/>		<b>Incr</b>	Min		Max			
Prozone Limit	0		<input type="radio"/>	Upper	<input type="radio"/>	Lower	Reagent ABS	<b>0</b>		<b>0.4</b>		
Unit	<b>mg/dL</b>		Decimal Point		<b>2</b>		Tech. Serum Lim.	<b>0</b>		<b>20</b>		
							Tech. Urine Lim.	<b>0</b>		<b>0</b>		
							Panic Limit	*		*		
Normal Values	AGE		Male		Female							
			Min	Max	Max	Max						
Serum	<b>Default</b>		<b>2.5</b>	<b>4.5</b>	<b>2.5</b>	<b>4.5</b>	Auto Dil.	<input checked="" type="radio"/>	<b>Yes</b>	<input type="radio"/>	No	
Serum							Y=aX+b	a =	<b>1</b>	b =	<b>0</b>	
Serum												
Urine values												
Calibration curve	<b>Straight</b>											

\* Data entered by the operator

In the **Std. Volume (Pre/Norm/Dil)** field, enter the **S. Vol. Normal ((Pre/Norm/Dil)** values