

# ALBUMIN

Spectrophotometric  
BROMOCRESOL GREEN

The application parameters comprised here constitute a guide to facilitate the validation of our reagents by the instrument. It is advisable to validate the use when there is any change in software or reagent versions.

## Instruments: EXPRESS 550 / Plus

### Reagent preparation

Reagent is ready to use.

### Instrument settings

Test Name	ALB	Test	ALB		
Test Bar Code					
Test Type	End Point	Curve Type	Blanked Linear		
Units	g/L	Nr. of Decimal Places	1		
Primary Wavelength	600	Secondary Wavelength	540		
Read Time/Interval	20	Sample Blank	No		
Factor					
Calibration Interval	168 hours	Nr. of Replicates	2		
Nr. of Calibrators	2				
Low Blank A Limit	-0.120	High Blank A Limit	1.000		
Low A Limit	-0.250	High A Limit	1.500		
Low Normal	35	High Normal	50		
Linearity Limit	70	Curve S.D. Limit	10.0		
Sample Volume	3 µL	Test	ALB		
Rerun Dilution Ratio	1				
Sample Diluent Bottle Type	*	Predilution Ratio	1		
Reagent Diluent Bottle Type	*				
	Reagent Volume	Bar Code	Diluent Volume	Lag Time	Bottle Type
Reagent 1	300 µL	AL1A		60	*
Reagent 2					

\* Enter the bottle type: plastic or glass